



(1) Resolution – swinging apertometer

Numerical aperture – measure of the ability of the lens to capture diffracted / scattered light.

Numerical aperture, $NA = n \sin \alpha$

n = refractive index of the medium between the object (specimen) and the objective lens

$n_{air} = 1.0003 \approx 1$



Angle 1	Angle 2	Collection angle	Half angle, α	Measured NA	Actual* NA	Resolution* (μm)	Depth of field* (μm)

- Calculate the N.A for: objective lens specimen

- What if we increase the collection angle to 164°?
- The Abbe formula to calculate lateral resolution is:

$$d = \lambda / 2NA$$

The approximate formula to calculate depth of field is:

$$\text{depth of field} = n\lambda / NA^2$$

- For calculations using white light, assume a wavelength of 500 nm (0.5 μm).